


# Structures System

## *Imagine • Design • Analyze*


### Experience real-world design

- ▶ Build bridges, cranes, hydraulic lifts, and more
- ▶ Structures can be reconfigured quickly; students can explore different designs
- ▶ Measure tension and compression



**Tied Arch Bridge**  
(See page 11)

**Load Cell & Amplifier Set**  
(See page 4)



The PASCO Capstone™ graph shows changes in the compression and tension in the supporting members as the car traverses the bridge.



**Materials Testing System**  
*Makes your study of structures complete.*

(See pages 14-15)

### Inside:

Measurement.....	Pages 2-3
Load Cells and Amplifiers.....	Page 4
Truss Set.....	Page 5
Bridge Set.....	Pages 6-7
Advanced Set .....	Pages 8-9
Human Structures .....	Page 10
Large Structures Set.....	Page 11
Classic Statics .....	Page 12
Hydraulics .....	Page 13
Materials Testing.....	Pages 14-15
Cast Beams .....	Page 16
Data Acquisition and Display .....	Page 17
Force and Displacement.....	Page 18
Bridges with Rigid Roadbeds .....	Page 19
Structures Resonance .....	Pages 20-21
Replacement Spares Sets .....	Pages 22-23

**PASCO**  
**2018**




# Structures System Overview

## Build bridges, cranes, catapults and roller coasters.

Experience real-world design building a large variety of structures. This reconfigurable system allows students to measure static and dynamic forces using load cells and still have time to redesign and test again.

Load Cell & Amplifier Set  
(See page 4)



As the car crosses the bridge, the forces measured by each load cell are graphed in real-time in PASCO Capstone™. Notice the diagonal member (green trace) switches from compression to tension as the car passes by.



**Only PASCO allows students to construct an endless array of structures and measure forces on them in real time.**



## Imagine... Design... Analyze...

Trusses, bridges, roller coasters, cranes, booms, human models and much more can all be quickly built and analyzed. Far more advanced than toothpick models and much more hands-on than computer simulations, PASCO Structures are ideal for real-world design.

And the ease of use of PASCO Structures Systems allows students to quickly build, test and then redesign their structures quickly and efficiently which supports the engineering process.



*Designing and building structures is simple and easy.*



*I-beams fit into connectors and are secured with thumb screws.*




*Load Cells measure the compression and tension...*



*...and may be placed anywhere in your structure.*



*Load Cells are then plugged into a Load Cell Amplifier which connects to a PASCO interface which is connected by USB or Bluetooth to a device running PASCO software.  
(See interfaces page 17)*



*PASCO software then allows students to view a graph of the forces in real-time and analyze the resulting data with a suite of powerful tools. PASCO Capstone runs on Mac or Windows Computers. SPARKvue software runs on computers, iPads, Android tablets, or Chromebooks.  
(See software page 17)*

# Load Cells and Amplifiers

## Choice of Load Cell Amplifiers:


### Load Cell Amplifier (6 ports)

PS-2198



This Load Cell Amplifier can accommodate up to six load cells and only needs a single PS-2100A USB Link (p. 2) to connect to a computer. Useful for doing an extensive analysis of a bridge by inserting six load cells at various points in the structure.

The Amplifier accepts either the 100N load cell or the 5N load cell or a combination of both. The maximum data sample rate is 500 Hz for each port.



Shown using the Truss Set ME-6960 and Load Cells.

**Load Cell Amplifier (6 ports)..... PS-2198**

#### Requires:

Load Cell 100N..... PS-2200  
Load Cell 5N..... PS-2201  
Interface..... (See page 17)

#### Also available:

Load Cell and Amplifier Set..... PS-2199

#### Set Includes:

Load Cell Amplifier  
Load Cell 100N (4)  
Instruction manual




### Dual Load Cell Amplifier

PS-2205



Use for applications where only one or two load cells are needed, such as measuring the force on the track at the top of a roller coaster loop. If you only want to examine the forces in a bridge one at a time, you can move a single load cell around in the bridge.

The Amplifier accepts either the 100N load cell or the 5N load cell or a combination of both. The maximum data sample rate is 1000 Hz for each port.



Shown using the Physics Set ME-6989.

**Dual Load Cell Amplifier..... PS-2205**

#### Requires:

Load Cell 100N..... PS-2200  
Load Cell 5N..... PS-2201  
Interface..... (See page 17)

#### Also available:

Load Cell and Dual Amplifier Set..... PS-2206

#### Set Includes:

Dual Load Cell Amplifier  
Load Cell 100N (1)




## Two Ranges of Load Cells:

### Load Cell 100N and 5N

Both load cells can be used with the same amplifier in any combination. The semi-transparent case lets students see the strain gauge and beam inside.

**Load Cell 100N..... PS-2200**  
**Load Cell 5N..... PS-2201**




Specifications		
	Load Cell 100N PS-2200	Load Cell 5N PS-2201
<b>Range:</b>	-100 N to +100 N	-5 N to +5 N
<b>Accuracy:</b>	±1% (± 1 N)	±1% (±0.05 N)
<b>Resolution:</b>	0.02 N	0.001 N
<b>Safe Overload:</b>	-150 N to +150 N	-7.5 N to +7.5 N

## Truss Set


ME-6990

- ▶ Teach the basics of trusses
- ▶ Demonstrate the properties of I-Beams


Use the Truss Set to build a variety of structures to investigate the principles of trusses. The ABS plastic I-Beams fasten securely together using the provided connectors and thumb screws. Load cells can be inserted anywhere into the design by replacing one beam at a time. Students can load the truss by hanging weights.



Students can construct a roof truss to study how the roof is supported in buildings.



Kingpost Roof Truss



Deck Truss Bridge

### Truss Set Includes:

One package each of Truss Set Members and Truss Set Screws.  
See pages 22-23 for details.

Truss Set ..... ME-6990

### Recommended:

Load Cell and Amplifier Set..... PS-2199  
includes four load cells. (See page 4))



# The Bridge Set

## Bridge Set


ME-6991

- ▶ Study the principles of bridge construction
- ▶ Build different scales
- ▶ Add Load Cells to see dynamic loading as car traverses bridge
- ▶ Design a roller coaster

Flexible road bed

Car with Mass

Warren with Verticals  
1.5 m long



The Bridge Set includes all the I-beams and connectors required to build the structures shown here. Special cord locks allow tensioning of cord (cables) for cross bracing. A flexible plastic road bed clips to the cross-beams and, using load cells, the tension and compression of each element can be displayed in real time.

## Build different types of bridges

Howe



Pratt




Warren



## Build different scale bridges


Students will learn the advantages of building taller bridges. Here, the larger bridge is constructed with a #4 I-beam on the diagonal, while the smaller bridge uses a #3 I-beam on the diagonal.




Both of these bridges can be built with the ME-6991 Bridge Set.

## Add load cells to measure static forces anywhere in the structure

Forces measured by Load Cells are displayed on a computer using PASCO Capstone™ Software (see page 17). A positive value represents compression.





I-Beams key into the load cell and are fastened with thumb screws.



## Design your own roller coaster

PASCO's Bridge Set allows students to design and build their own roller coaster for detailed studies of conservation of energy and centripetal force. The flexible track is perfect for building hills, valleys and even a loop! Car with low-friction ball-bearing wheels minimizes energy losses. Measure the speed of the car using photogates or a Motion Sensor.



### Bridge Set Includes:

Two packages each of Truss Set Members and Truss Set Screws  
One package each of Roadbed Spares and Cord Lock Spares

See pages 22-23 for details.

**Bridge Set**..... **ME-6991**

### Recommended:

Load Cell and Amplifier Set..... PS-2199  
includes four load cells. (See page 4))

### Shown in use with:

Accessory Photogate..... ME-9204B

### Requires:

Interface..... (See page 17)








# The Advanced Structures Set

## Advanced Structures Set

ME-6992B

- ▶ Build larger bridges
- ▶ Build cranes, catapults, cars and


The Advanced Structures Set includes more components to build larger structures. Axles and pulleys allow construction of cranes, cars and even a working catapult!






## Human Leg Model

The articulated leg shown below uses a rubber band (not included) for the quadriceps and has a load cell on the foot to measure the force that the “toe” exerts on the ball. The impulse (area under the curve) is equal to the resulting momentum of the ball.





See page 10 for more Human Structures.



## Rubber Band Car


Build a working catapult, a car with rubber band suspension, and a rubber band-powered car. The “rubber band” car shown here uses the ME-8986 Rubber Cord (sold separately) and travels over 50 ft. in 10 seconds.



Use these flexible I-Beams to make a bridge which dramatically demonstrates how a bridge fails and yet the beams will return to their original shape when the load is removed.



For more examples using the Advanced Structures Set, see pages 8, 17, 18 and 19.



### Advanced Structures Set Includes

Truss Set Members (3 pkgs.)	Round Connector Spares (1 pkg.)
Truss Set Screws (4 pkgs.)	Angle Connector Spares (1 pkg.)
Flexible I-Beams (1 pkg.)	Flat Beams (1 pkg.)
Cord Lock Spares (1 pkg.)	Structures Rod Clamps (1 pkg.)
Axle Spares (1 pkg.)	Force Platform Structures Bracket

See pages 22-23 for details.

**Advanced Structures Set**.....**ME-6992B**

### Shown in use with:

Load Cell and Amplifier Set  
includes four load cells. (See page 4)).....PS-2199

Additional Load Cell 100N.....PS-2200

Hooked Mass Set.....SE-8759

Large Slotted Mass Set.....ME-7566

Motion Sensor.....PS-2103A

Rubber Cord (spool of 30 m).....ME-8986

### Requires:

Interface.....(See page 17)

# The Human Structures Set

## Human Structures Set


ME-7001

- Build models that represent real life examples.
- Bring homework problems to life.

### Human Arm Model

Build a realistic arm model and directly measure the forces exerted by the biceps muscle (tension in supporting cord).


Vary the length and angle of upper and lower arm, as well as the point of attachment of the muscle.



Support Structure allows the angle of the upper arm to be easily adjusted.

Load cell represents the biceps muscle.

Axle allows arm to pivot freely.




Construct all three models concurrently with this set.

### Human Back Model

Model the forces acting on a human back. Vary all parameters including position of back muscle attachment and angle of the torso. Directly measure the force exerted by the back muscles.

Lumbar Muscles

Load cells directly measure forces exerted on back model.




#### Human Structures Set Includes

Five packages of Truss Set Screws  
Two packages of Truss Set Members  
Two packages of Connector Spares  
One package each of #6 I-Beam Spares, Cord Lock Spares, Axle Spares, Round Connector Spares, Angle Connector Spares, roll of rubber cord.


See pages 22-23 for details.

### Human Leg Model

Measure the force needed to support the leg at various angles.



Load cell represents quadriceps muscles.



Directly measure the force needed to support the leg at various angles.

Human Structures Set..... ME-7001

#### Shown in use with:

Load Cell and Amplifier Set..... PS-2199  
(includes four load cells)

Hooked Mass Set..... SE-8759


Large Slotted Mass Set..... ME-7566




## Large Structures Set

ME-7003

The Large Structures Set includes all the components contained in the Advanced Structures Set (ME-6992B) plus additional parts to build even bigger structures. It also includes the Mini Cars with plastic track to build roller coasters and to add realistic roadbeds to your bridges.



**Double Tied Arch Bridge**  
2.8 m long



**Arch Truss**  
2 m long



**Cable Stayed**  
3.8 m long




**Suspension Bridge**  
3 m long



*Build this crane and 15 other HUGE structures with this one comprehensive set!*

*Add load cells (page 4) to measure forces anywhere in the structure.*

*Student uses the Hydraulic/Pneumatic System (page 13) to adjust the crane.*



**Drawbridge**  
3.2 m long

*Student uses the Hydraulic/Pneumatic System (page 13) to raise and lower the drawbridge.*

### Large Structures Set Includes

Six packages of Truss Set Screws  
Three packages of Truss Set Members  
Two packages of Connector Spares  
One package each of #6 I-Beam Spares, Flexible I-Beams, Cord Lock Spares, Axle Spares, Round Connector Spares, Angle Connector Spares, Flat Beams, Structures Rod Clamps, Mini Car Track Spares, Force Platform Structures Bracket, and one each Green Car, Yellow Car, 9.1 m Track, and Starter Bracket

**Large Structures Set** ..... ME-7003

### Shown in use with:

Load Cell and Amplifier Set  
includes four load cells..... PS-2199  
Hydraulic/Pneumatic Structures ..... ME-6984  
Slotted Mass Set..... ME-7589


### Requires:

Interface..... (See page 17)


# Classic Statics Using the Advanced System

## Forces on a Boom


Vary all parameters including length and angle of the boom. Directly measure the horizontal and vertical forces exerted by the pivot (axle) on the boom, and the tension in the supporting cord.



Support structure is built using parts from the Advanced Structures Set, and uses 1/2 kg masses from the ME-7566 Large Slotted Mass Set (sold separately) for counter balance.




Classic statics problem of pivoted boom supported by a cable.



## Teeter Totter

Take "meter stick" torque to a new level! By building their own unique structures, students learn about center of mass, torque, and static equilibrium as never before.




## A Lesson in Balance

Circus performer not included!




## Ladder Against Wall

How does the required frictional force (measured by the load cell at bottom) change as the location of the hanging mass is moved further up the ladder?



## Stability and Center of Mass



The structures shown here are built with a foam core center (not included). Each object is hung from several different pivot points to locate the center of mass. Is the structure on the right stable or unstable?

Advanced Structures Set..... ME-6992B

**Shown in use with:**

Load Cell and Amplifier Set..... PS-2199

(includes four load cells)

Additional 100 N Load Cell ..... PS-2200

Hooked Mass Set..... SE-8759

Large Slotted Mass Set..... ME-7566

Angle Indicator ..... ME-9495A

Mass and Hanger Set..... ME-8979

Large Table Clamp ..... ME-9472

Steel Rod (90 cm)..... ME-8738

**Requires:**


Interface..... (See page 17)



## Hydraulic and Pneumatic Structures

ME-6984


Add a hydraulic/pneumatic ram to make your structures move and do work. Not only will students see the cranes and jacks in action, they can directly measure the pressure and volume to calculate how much work was done.



An Ideal Gas Law Apparatus (TD-8596A), which has an internal thermistor, is used to pump air into the cylinder. A Pressure/Temperature Sensor (PS-2146) records the air pressure and temperature while the Rotary Motion Sensor (PS-2120) records the movement.



Valves are used with the syringe to pump up this fork lift. The use of different size syringes shows how a smaller force requires a greater number of pumps to do the same amount of work as a larger force.



The weight is lifted using a syringe of water to fill the master cylinder. An Absolute Pressure Sensor (PS-2107) measures the pressure and a Rotary Motion Sensor (PS-2120) records the movement of the piston.



This scissor lift uses pulleys to change the mechanical advantage.



**Includes**  
 Master Cylinder  
 Pressure Sensor "T"  
 Check Valves and Tubing  
 10 ml Syringe  
 20 ml Syringe  
 60 ml Syringe  
 Drive belt for Rotary Motion Sensor (not shown)

### Hydraulic/Pneumatic Structures..... ME-6984

#### Required:

Advanced Structures Set.....	ME-6992B
Steel Rod (90 cm).....	ME-8738
Absolute Pressure Sensor.....	PS-2107
Rotary Motion Sensor.....	PS-2120A
Pressure/Temperature Sensor.....	PS-2146
Ideal Gas Law Apparatus.....	TD-8596A

#### Not shown but required for data collection:

Interface and PASCO Capstone™ Software ..... (see page 17)

# Materials Testing


## Materials Testing Machine

ME-8236


- ▶ 7100 N max load
- ▶ Hand-cranked so students can feel samples break
- ▶ Inexpensive samples make it possible for each student to experience it firsthand

Measure force and displacement for various materials as they are stretched, compressed, sheared, or bent. Investigate material properties including Young's Modulus, Tensile Strength, Yield Strength, Ductility and Modulus of Resilience.

The Materials Testing Machine measures force with a 7100 N load cell and displacement with an optical encoder. It runs on PASCO Capstone software, which has a built-in compliance calibration wizard and has all the tools to record and display stress vs. strain, apply linear fits to find Young's Modulus, and to record and play back webcam movies of the breaking samples synced to the data.



Tensile stress versus strain is plotted in PASCO Capstone software for steel, annealed steel, and brass.



For annealed steel, a linear fit is applied to find Young's Modulus.

### Specifications

Load cell capacity: 7100 N (1600 lbs)

Machine weight: 20 lbs (9 kg)

Footprint: 24 wide x 25 depth x 51 cm height

Lead screw length: 38 cm

Sturdy base: cast aluminum

Mounting holes: for bolting to table

### Materials Testing Machine Includes


Machine

Compliance Calibration Rod

Safety Shields

Requires Capstone software (See page 17)

(does not include any samples)



Shown in use without the included safety shields.



Tensile Samples (set of 10 each)  
Shown using tensile sample from Steel Set (ME-8243).

**Materials Testing Machine**..... **ME-8236**

**Tensile Samples** (set of 10 each)

Aluminum..... ME-8231

Brass..... ME-8232

Annealed Steel..... ME-8233

Steel..... ME-8243

Acrylic..... ME-8234

Polyethylene..... ME-8235

### Required:

PASCO Capstone™ Software ..... (See page 17)



## Comprehensive Materials Testing System

ME-8244

- ▶ Compression and tensile testing
- ▶ Column buckling
- ▶ Three-and four-point bending
- ▶ Shear testing
- ▶ Stress lines with photoelasticity


System includes everything needed to study material testing: Testing Machine with attachments, test samples, and computer software. PASCO Capstone Workbooks include setup instructions, theory, and detailed analysis questions.



Perform three-point and four-point bending.



Perform a four-point bend test on the ME-6983 Cast Beams from the PASCO Structures System. Quantities measured include the Flexural Elastic Modulus and the Modulus of Rupture for the material.



Shown performing tensile test of brass sample.



Visualize the stress lines as the clear photoelastic beam is compressed and buckles.



### Comprehensive Materials Testing System Includes

Materials Testing Machine

Storage Base

Tensile Samples (10 of each): aluminum, brass, annealed steel, steel, acrylic, and polyethylene

Bending Accessory

Four-point Bending Load Anvil

Photoelasticity Accessory (with photoelastic beams)

Shear Accessory (with Shear Samples)

Structures Beam Fixture

Thin Beams

Cast Spares

Compression Accessory (with Compression Samples)

Flat Coupon Fixture

Plastic Flat Coupons

Metal Flat Coupons

Clevis Grip

10-32 Adapter

AirLink Interface

PASCO Capstone Software Single User License (See page 17)




Comprehensive Materials Testing System ..... ME-8244

# Cast Beams Set

## Cast Beam Structures Set


ME-7009

Make your own cast beams that look like pre-stressed concrete beams. Test them and you'll find they perform like them, too. These beams are cast with a mixture of sand and plaster of Paris (not included). The rebar is made of the same plastic used for the I-beams. Students can explore how the strength of the beam is affected by the amount of tension put on the rebar, the mixture of sand and plaster of Paris, or using one or two rebar.




Both the tension fixture and the test fixture can be built concurrently with this set.

The graph of hanging mass versus displacement shows the relative strengths of three beams: one cast beam made with no pre-load; one cast beam made with 60N of pre-load; and one normal plastic I-beam. Notice that the traces for the cast beams show discontinuities when the beams cracked. Also notice that the pre-loaded cast beam is stronger than the plastic I-beam until the cast beam cracks.




Tension Fixture




Displacement of beam is measured as load is increased.

The beam rebar is kept under tension while the beam is drying.



**Step 1:** The rebar with connecting ends snaps into the plastic mold. Pour a mixture of sand and plaster of Paris into the mold.




**Step 2:** After it dries, it is easy to remove the flexible plastic mold from the cast beam.

The cast beam shown here is tested to destruction under compression using the Materials Testing Machine on pages 14-15.



Close-up of beam after destruction showing plastic "rebar."



PASCO Capstone graph shows that strength of Plaster of Paris increases with cure time.

## Cast Beam Spares

Consumable replacement parts for Cast Beams; these can also be used with the Advanced Structures Set (page 8).



### Includes

10 Reusable Plastic Molds  
30 Rebar with Connectors

Cast Beam Spares ..... ME-6983

### Cast Beam Structures Set Includes

One package each of Truss Set Members, Cord Lock Spares, Axle Spares, Round Connector Spares, Angle Connector Spares and two packages of Truss Set Screws  
See pages 22-23 for details.

Cast Beam Structures Set .....ME-7009

### Also shown:

Displacement Sensor .....PS-2204  
Large Slotted Mass Set .....ME-7566  
Round Base with Rod .....ME-8270

### Not shown but required for data collection:

Interface and  
PASCO Capstone™ Software .....(see next page)



## PASCO's Data Collection and Analysis Software

### PASCO Capstone™ Software

- Compatible with all PASCO USB interfaces
- For use with Mac® and Windows® computers

Whether you want your students to explore and create lab write-ups on their own, or you want to tailor a lab write-up with very specific instructions, PASCO Capstone™ has the power and flexibility to meet the needs of your lab.



**PASCO Capstone Single User License**  
for Mac®/Windows®..... **UI-5401**  
**PASCO Capstone K-12 Campus License**  
for Mac®/Windows®..... **UI-5405**

For more information visit [www.pasco.com/capstone](http://www.pasco.com/capstone)

### SPARKvue® Software

- For use with mobile devices or Chromebooks
- Free app for iPads, Android tablets and Chromebooks



**SPARKvue**  
(single user license)..... **PS-2401**  
(site license) ..... **PS-2400**

Free mobile apps. Visit [pasco.com/downloads](http://pasco.com/downloads)

For more information visit [www.pasco.com/sparkvue](http://www.pasco.com/sparkvue)

## Connect PASPORT Sensors to a Computer



AirLink, has one PASPORT sensor port and is compatible with our full line of sensors. It has USB and Bluetooth™ 4.0 connectivity and comes with a USB cable.



**AirLink** ..... **PS-3200**

## Fast, Flexible and Powerful!

### The PASCO 550 Universal Interface

- 1 MHz sampling rate
- 2 high-speed analog inputs
- 2 digital inputs for photogates and other timing sensors
- 2 PASCO PASPORT sensor inputs
- Signal generator with built-in Voltage and Current sensors
- Runs on Capstone or SPARKvue
- Bluetooth® connectivity



**550 Universal Interface**..... **UI-5001**

**Requires:**


PASCO Capstone Software **OR** SPARKvue Software

## The Ultimate Sensor Interface for Physics and Engineering

### 850 Universal Interface

- Rugged design
- Fully compatible and expandable
- An incredible value
- Runs on PASCO Capstone™

Here's the most powerful educational lab interface in the world. Compatible with over 120 PASCO PASPORT® and ScienceWorkshop® Sensors.



**850 Universal Interface**..... **UI-5000**

**Requires:**

PASCO Capstone Software **OR** SPARKvue Software

For more information visit [www.pasco.com/850](http://www.pasco.com/850)


# Force and Displacement

## Measure support forces with a Force Platform

### Force Platform

PS-2141

Measure the support forces of a crane by connecting it to a Force Platform (PS-2141) using the special Force Platform Structures Bracket (ME-6988). The Force Platform is supported by four individual load cells that combine to measure the total vertical force on the platform. These four readings can also be viewed separately to measure the unequal forces on the crane supports.




Crane built using the Advanced Structures Set (ME-6992A) shown on page 8.

**Force Platform.....PS-2141**  
**Force Platform Structures Bracket.....ME-6988A**

## Measure passive damping with a Motion Sensor

This building frame is built with an Advanced Structures Set using the Flat Beams. A pendulum with drag caused by strings is suspended from the top of the building. The Motion Sensor is positioned to record the oscillation of the building.



The gray graph, in PASCO Capstone™ Software (page 17) shows the oscillation without the pendulum. The red graph shows the damping caused when the pendulum is allowed to oscillate.

**Advanced Structures Set.....ME-6992B**  
**Motion Sensor.....PS-2103A**  
**Large Slotted Mass Set.....ME-7566**

## Measure bridge deflection with a Displacement Sensor

### Displacement Sensor

PS-2204

The Displacement Sensor measures the travel of a spring-loaded indicator pressed against a bridge as the bridge is loaded. It consists of a PASPORT sensor that plugs into the included Digital Indicator and a digital travel indicator that has its own digital LED readout and can be used as a stand-alone device. When the PASPORT sensor is plugged into an interface, the reading can be recorded.



### Specifications

Maximum Travel: 10 mm  
Maximum Sample Rate: 5 Hz  
Resolution: 0.013 mm (0.0005 in)



**Displacement Sensor.....PA-2204**

### Shown in use with:

Hooked Mass Set.....SE-8759  
Small "A" Base.....ME-8976  
60 cm long Steel Rod (threaded).....ME-8977

### Required:

Interface.....(See page 17)



**Combine the rigid, plastic PASTrack sections with components from the Structures System to build truss and cable-stayed bridges with realistic rigid decks.**

## PASTrack Cable-Stayed Bridge


The cable-stayed bridge shown here was constructed using components from the Advanced Structures Set and two ME-6997 Round Connector Spares sets. The roadbed uses four sets of PASTrack and four sets of curved PASTrack. The bridge is built starting with the center column and working symmetrically outward, always keeping the bridge in balance. As each 50 cm section of rigid deck is added to each side, new supporting cable is added.

**Cable-Stayed Bridge**  
1.5 m Tall

PASCO Capstone™ graph shows the change in tension and compression of the supporting members as the Motorized Cart moves across the span.



Add Load Cells to measure the support forces in the bridge.



## PASTrack Truss Bridge


Combine the plastic PASTrack sections with the components from the Advanced Structures Set (page 8) to build a rigid deck to support dynamics carts. Use load cells (page 4) to directly measure the forces as the cart traverses the bridge.

Adjustable Endstop

Components from Advanced Structures Set

PASTrack



Load Cell & Amplifier Set



Load cells are used to measure the tension and compression in the members.

## Dynamics Track Support Reactions

Combine the plastic PASTrack sections with the components from the Advanced Structures Set (page 8) to measure the support reactions as the Motorized Cart climbs the hill. The Load Cell and Amplifier Set (page 4) directly measures the forces in the structure as the Motion Sensor measures the cart's position.



### Equipment Shown:

**Advanced Structures Set ..... ME-6992B**  
 PAScar (set of 2) ..... ME-6950  
 PASTrack (2 sections) ..... ME-6960  
 Curved PASTrack ..... ME-6841  
 Motorized Cart ..... ME-9781  
 PAScar Cart Mass (set of 2) ..... ME-6757A  
 Adjustable Endstop ..... ME-8971  
 Round Connector Spares ..... ME-6997  
 Load Cell & Amplifier Set  
 (includes four load cells) ..... PS-2199  
 Motion Sensor ..... PS-2103A  
 PASCO Capstone Software ..... (see page 17)  
 Interface ..... (see page 17)



## Structures Resonance

PASCO's Structures System is perfect for demonstrating resonance in complex systems. The plastic I-Beams clearly show two different bending moments and can be connected together to build a variety of structures.




Composite  
I-Beam  
1.2 m tall

The long plastic I-Beam is constructed of components from the Advanced Structures Set (page 8). It is driven using the Mechanical Wave Driver and the Function Generator, demonstrating the three lowest harmonics.

### Equipment Shown:


**Advanced Structures Set** ..... ME-6992B  
Function Generator ..... PI-8127  
Mechanical Wave Driver ..... SF-9324  
Large Slotted Mass Set ..... ME-7566  
5N Load Cell ..... PS-2201  
45cm Stainless Steel Rod ..... ME-8736  
Large Rod Base ..... ME-8735

Shaking Tower  
75 cm tall



Mechanical  
Wave Driver  
(SF-9324)

See page 18 for more  
examples using the  
Shaking Tower.




Function  
Generator  
(PI-8127)



## Measure acceleration with a 5N Load Cell


Connect one end of a load cell to the structure and attach a mass to the other end of the load cell. The acceleration of the structure is measured in real time as the structure shakes.



## Bridge Vibrations

Concepts:

- ▶ Study resonance in complex systems
- ▶ Compare driven vs. free vibrations



5N Load Cells are used to measure the oscillations of the bridge at different positions.

The bridge is struck by hand and allowed to freely oscillate. The FFT (using PASCO Capstone™) shows several resonant frequencies.

The resonance of the bridge is characterized by driving the bridge at different resonant frequencies. Note how different the amplitudes are at different locations on the bridge.

### Download This Experiment

The FREE experiment files include instructions in Microsoft Word™, PASCO Capstone™ workbook files with sample data, and graphics. Download these experiments at [www.pasco.com/CapstoneExperiments](http://www.pasco.com/CapstoneExperiments).

#### Experiment Includes

Large Structures Set  
Load Cell Amplifier  
100N Load Cell  
5N Load Cell (5)  
Mechanical Wave Driver

ME-7003  
PS-2198  
PS-2200  
PS-2201  
SF-9324

4 mm Banana Plug Cords  
Rubber Cord  
Large Slotted Mass Set (4)  
Short Mass Hanger (2)  
20 g Masses (3 sets of 6)

SE-9750  
ME-8986  
ME-7589  
ME-7590  
ME-8983

#### Bridge Vibrations ..... EX-5548

##### Required:

850 Universal Interface (see page 17)  
PASCO Capstone Software (see page 17)

# Replacement Spares Sets

## Truss Set Members

### Includes:

I-beam #5 (8) 24 cm long  
I-beam #4 (18) 17 cm long  
I-beam #3 (18) 11.5 cm long  
I-beam #2 (8) 8 cm long  
I-beam #1 (8) 5.5 cm long  
Connectors (14)

**Truss Set Members.....ME-6993**

## Truss Set Screws

### Includes:

75 screws  
All components in the Structures System use this same 6-32 thumb screw.

**Truss Set Screws .....ME-6994**

## Connectors Spares

Set of 14 connectors used to join truss members. This is the same connector included in the Truss Set ME-6990.

**Connector Spares .....ME-7002**

## Angle Connector Spares

### Includes:

Sliding connectors (12),  
Angle connectors (24),  
Straight connectors (24).

### Angle Connector

**Spares.....ME-6999A**

## Round Connector Spares

### Includes:

Round connectors (6)  
Flat connectors (6)  
Six bolts with nuts.

### Round Connector

**Spares.....ME-6997**

## Cord Lock Spares

### Includes:

32 cord-tensioning clips and a spool of yellow cord.

**Cord Lock Spares.....ME-6996**

**Yellow Cord (2 pack).....ME-9876**

## Axle Spares

### Includes:

Drive wheel with rubber tire (4), pulleys with "O" rings (12 each), axles (two each of three lengths), spacers (12) and collets (24).

**Axle Spares .....ME-6998A**

## Thin I-Beams

### Includes:

Thin I-beam #4 (24) 17 cm long  
Thin I-beam #3 (24) 11.5 cm long

**Thin I-Beams .....ME-7012**

## Photoelastic Beams

### Includes:

Clear, Polycarbonate Thin I-beams #4 (24) 17 cm long, and #3 (24) 11.5 cm long.

**Photoelastic Beams.....ME-7011**

## Flexible I-Beams

### Includes:

Flex I-beam #5 (10) 24 cm long  
Flex I-beam #4 (18) 17 cm long  
Flex I-beam #3 (18) 11.5 cm long

**Flexible I-Beams .....ME-6985**

## Flat Beams

### Includes:

2x3 beams 12 cm long  
F4 beams 17 cm long  
3x4 beams 19 cm long

**Flat Beams .....ME-6987**

## #6 I-Beam Spares

Longer beam to supplement the Truss Set ME-6990.

Has the same cross section as the shorter beams.

### Includes:

24 of the #6 I-beams, 35 cm long.

**#6 I-Beam Spares.....ME-7008**

## Cast Beam Spares

Consumable replacement parts for Cast Beams ME-7009. Includes 10 Reusable Plastic Molds and 30 Rebar with Connectors. These can also be used with the Advanced Structures Set ME-6992B.

**Cast Beam Spares .....ME-6983**



## Roller Coaster Track (9.1 m)



Longer replacement roll of flexible plastic track for use with the Bridge Set ME-6991, Physics Structures Set ME-6989, Large Structures Set ME-7003, and Roadbed Spares ME-6995.

**Roller Coaster Track .....ME-9814**

## Roadbed Spares

Starter Bracket



### Includes:

Flexible roadbed (3 m)  
Roadbed clips (24)  
Car with flag  
Extra mass, mini car starting bracket, and track couples (2).

**Roadbed Spares.....ME-6995**

## Force Platform Structures Bracket



### Includes:

Brackets (2)  
Screws (4)

**Force Platform  
Structures Bracket .....ME-6988A**

## Mini Cars (Set of 3)

These cars feature low-friction ball bearings and ABS construction to withstand repeated impacts. One red, one yellow, and one green Mini Car included. Each car includes a slot for a supplied photogate flag, cup/mass holder, and cup. The body of the car extends just far enough below the wheels to protect them should the car leave the track.



### Includes:

Mini Cars (3)  
Decals  
Elastic Bumpers (3)  
Flags (3)  
Rubber Bands

**Mini Cars (Set of 3) .....ME-9813**

## Mini Car Track Spares



### Includes:

Two gates  
Two track couplers  
One bag (24) of roadbed clips

**Mini Car Track Spares .....ME-6974**

## Large Slotted Mass Set

Consists of nine iron disks of 0.5 kg each. The hanger is 0.5 kg. Each piece cast and machined to 1 gram accuracy. The 2 kg Mass Set is similar in mass and hanger weight to the Large Slotted Mass Set, but this more compact set includes a shorter 1/2 kg hanger and three 1/2 kg slotted masses. The size and weight of these sets make both ideal for creating loads on the Structures Systems.



**Large Slotted  
Mass Set.....ME-7566**

## Hydraulic and Pneumatic Structures



### Includes:

Master Cylinder  
Pressure Sensor "T"  
Check Valves and Tubing  
Syringes (10, 20, 60 ml)  
Drive Belt for Rotary Motion Sensor (not shown)

**Hydraulic/Pneumatic  
Structures.....ME-6984**

## Structures Rod Clamps (Set of 2)

Connects  
structure  
members to  
1/2 inch rod.



**Structures  
Rod Clamps (2) .....ME-6986**

# PASCO

10101 Foothills Blvd. • Roseville, CA 95747-7100

**+1 916-462-8383**

**[www.pasco.com/structures](http://www.pasco.com/structures)**



## Dystrybutor PASCO w Polsce:

**image**

IMAGE RECORDING SOLUTIONS Sp. z o.o.  
ul. Arkuszowa 190, 01-934 Warszawa  
tel. +48 22 752 27 88 – 96, fax. +48 22 752-27-97  
NIP:526-10-25-668, REGON:010994837

**PASCO**

**Zapraszamy na strony:** [www.pasco.com.pl](http://www.pasco.com.pl), [www.pasco.com](http://www.pasco.com)

**Prosimy o kontakt:** E [pasco@irs.com.pl](mailto:pasco@irs.com.pl), M +48 606 850 155, T +48 22 752 27 88 – 96 wew. 109

Profesjonalną i odpowiedzialną współpracą z nauczycielami wspieramy rozwój uczniów i studentów w współczesnym świecie.

## *Connect your students to the real world!*



This model of the Louisville 2nd Street Bridge was built using the PASCO Structures System. The section shown here is approximately 4 m long.